

Human TADA3L protein, His tag

Cat. No. GTX68158-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 49 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 200mM NaCl, 50% Glycerol, 2mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMSELKDC PLQFHDFKSV DHLKVCPRYT
AVLARSEDDG IGIEELDTLQ LELETLSSA SRRLRVLEAE TQILTDWQDK KGDRRFLKL RDHELGAPPK HGKPKKQKLE
GKAGHGPGPG PGRPKSKNLQ PKIQEYFTD DPIDVPRIK NDAPNRFWAS VEPYCADITS EEVRTLEELL KPPEDAEHY
KIPPLGKHYS QRWAQEDLLE EQKDGAARAA VADKKKGLMG PLTELDTKDV DALLKKSEAQ HEQPEDGCPF GALTQRLQA
LVEENIISPM EDSPIDMSG KESGADGAST SPRNQNKPFV VPHTKSLESR IKEELIAQGL LESEDRPAED SEDEVLAELR KRQAEKALS
AHNRTKKHDL LRLAKEEVSQ QELRQVRMA DNEVMDAFRK IMAARQKKRT PTKKEKDQAW KTLKERESIL KLLDGL

Expression System E. coli**Purity** > 80% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68158-pro Image

3μg Human TADA3L protein (GTX68158-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).