

## Human POLR3F protein, His tag

Cat. No. GTX68183-pro

|         |       |
|---------|-------|
| Species | Human |
|---------|-------|

|         |        |
|---------|--------|
| Package | 100 µg |
|---------|--------|

## Applications

|                   |         |
|-------------------|---------|
| Observed MW (kDa) | 36 kDa. |
|-------------------|---------|

## Properties

|      |        |
|------|--------|
| Form | Liquid |
|------|--------|

|        |                                                  |
|--------|--------------------------------------------------|
| Buffer | 20mM Tris-HCl, 200mM NaCl, 30% Glycerol, 1mM DTT |
|--------|--------------------------------------------------|

|              |                  |
|--------------|------------------|
| Preservative | No preservatives |
|--------------|------------------|

|         |                                                                                                                                                                                                                            |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|               |                                                                            |
|---------------|----------------------------------------------------------------------------|
| Concentration | 0.5 mg/ml (Please refer to the vial label for the specific concentration.) |
|---------------|----------------------------------------------------------------------------|

## Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAEVKVK VQPPDADPVE IENRIIELCH QFPHGITDQV IQNEMPHIEA QQRAVAINRL LSMGQLDLLR SNTGLLYRIK DSQNAGKMKG SDNQEKLVYQ IIEDAGNKGI WSRDIRYKSN LPLTEINKIL KNLESKKLIK AVKSVAASKK KVMYLYNLQP DRSVTGGAWY SDQDFESEFV EVLNQQCFKF LQSKAETARE SKQNPMIQRN SSFASSHEVW KYICELGISK VELSMEDIET ILNTLIYDGK VEMTIIAAKE GTVGSDVGHM KLYRAVNPII PPTGLVRAPC GLCPVFDDCH EGGEISPSNC IYMTIEWLEF

|                   |         |
|-------------------|---------|
| Expression System | E. coli |
|-------------------|---------|

|        |                    |
|--------|--------------------|
| Purity | > 90% by SDS-PAGE. |
|--------|--------------------|

|             |              |
|-------------|--------------|
| Conjugation | Unconjugated |
|-------------|--------------|

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX68183-pro Image

3μg Human POLR3F protein (GTX68183-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).