

Human HSP105 protein, His tag

Cat. No. GTX68194-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 97 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 50mM NaCl**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMSV GLDVGSQSCY IAVARAGGIE TIANEFSDRC TPSVISFGSK NRTIGVAAKN QQITHANNTV SNFKRFHGRA FNDPFIQKEK ENLSYDLVPL KNGGVGKVM YMGEELFSV EQITAMLLTK LKETAENSLK KPVTDCVISV PSFTDAERR SVLDAAQIVG LNCLRLMNDM TAVALNYGIY KQDLPSLDEK PRIWVVDMDG HSAFQVSACA FNKGKLVLG TAFDPFLGGK NFDEKLVEHF CAEFKTKYKL DAKSKIRALL RLYQECEK KLMSSNSTDL PLNIECFMND KDVSGKMNRS QFEELCAELL QKIEVPLYSL LEQTHLKVED VSAVEIVGGA TRIPAVKERI AKFFGKDIST TLNADEAVAR GCALQCAILS PAFKVRFSV TDAVPFISL IWNHDSDETE GVHEVFSRNH AAPFSKVLTF LRRGPFELEA FYSDPQGVYV PEAKIGRFV QNVSAQKDGE KSRVKVKVRV NTHGIFTIST ASMVEKVPTE ENEMSSEADM ECLNQRPPEN PDTDKNVQQD NSEAGTQPQV QTDAAQTSQS PPSPELTSEE NKIPDADKAN EKKVDQPPEA KKPPIKVVNV ELPIEANLVW QLGKDLLNMY IETEGKMIMQ DKLEKERND KNAVEEYVYE FRDKLCGPYE KFICEQDHQN FLRLLTETED WLYEEGEDQA KQAYVDKLEE LMKIGTPVKV RFQEAERPK MFEELGQRLQ HYAKIAADFR NKDEKYNHID ESEMKKVEKS VNEVMEWMNN VMNAQAKKSL DQDPVWRAQE IKTKIKELNN TCEPVVTQPK PKIESPKLER TPNGPNIDKK EEDLEDKNNF GAEPHQNGE CYPNEKNSVN MDLD

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 09 Page 1 of 2

DATA IMAGES



GTX68194-pro Image

3μg Human HSP105 protein (GTX68194-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).