

Human SPIN1 protein, His tag

Cat. No. GTX68203-pro

Species	Human
---------	-------

Package	20 µg
---------	-------

Applications

Observed MW (kDa)	30 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 150mM NaCl, 20% Glycerol
--------	---

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	---

Region/Sequence	Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMKTPFGK TPGQSRADA GHAGVSANMM KKRTSHKKHR SSVGPSKPVS QPRRNIVGCR IQHGWKEGNG PVTQWKGTVL DQVPVNPSLY LIKYDGFDCV YGLELNKDER VSALEVLPRD VATSRISDAH LADTMIGKAV EHMFTEDGS KDEWRGMVLA RAPVMNTWYF IYKDPVLY MYQLDDYKE GDLRIMPDSN DSPPAEREPG EVVDSL VGKQ VEYAKEDGSK RTGMVIHQVE AKPSVYFIKF DDDFHIVYVD LVKTS
-----------------	--

Expression System	E. coli
-------------------	---------

Purity	> 90% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------


Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---

Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68203-pro Image

3μg Human SPIN1 protein (GTX68203-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).