

Human MRG15 protein, His tag

Cat. No. GTX68204-pro
Species Human

Package
100 µg

Applications
Observed MW (kDa) 41 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 150mM NaCl, 30% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSHMAPKQD PKPKFQEGER VLCFHGPLY
EAKCVKVAIK DKQVKYFIHY SGWNKNWDEW VPESRVLYV DTNLQKQREL QKANQEYAE GKMRGAAPGK KTSGLQQKNV
EVKTKKNKQK TPGNGDGGST SETPQPPRKK RARVDPTVEN EETFMNRVEV KVKIPEELKP WLVDWDLIT RQKQLFYLP
KKNVDSILED YANYKKS RGN TDNKEYAVNE VVAGIKEYFN VMLGTQLLYK FERPYAEIL ADHPDAPMSQ VYGAPHLLRL
FVRIGAMLAY TPLDEKSLAL LLNYLHDFLK YLAKNSATLF SASDYEVAPP EYHRKAV

Expression System E. coli

Purity > 85% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68204-pro Image

3 μ g Human MRG15 protein (GTX68204-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).