

Human TOMM34 protein, His tag

Cat. No. GTX68206-pro

Species	Human
---------	-------

Package	100 µg
---------	--------

Applications

Observed MW (kDa)	35 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM DTT
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	---

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAPKFPD SVEELRAAGN ESFRNGQYAE
ASALYGRALR VLQAQGSSDP EEESVLYSNR AACHLKDGNCRDCIKDCTSA LALVPFSIKP LLRRASAYEA LEKYP MAYVD
YKTVLQIDDN VTSAVEGINR MTRALMDSL GPEWRLKLPSI PLVPVSAQKR WNSLPSENHK EMAKSKSKET TATKNRVPSA
GDVEKARVLK EEGNELVKKG NHKKAIEKYS ELLCSNLES ATYSNRALCY LVLKQYTEAV KDCTEALKLD GKNVKAFYRR
AQAHKALKDY KSSFADISNL LQIEPRNGPA QKLRQEVKQN LH

Expression System	E. coli
-------------------	---------

Purity	> 85% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68206-pro Image

3μg Human TOMM34 protein (GTX68206-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).