

Human EXOSC8 protein, His tag

Cat. No. GTX68251-pro
Species Human

Package

50 µg

Applications
Observed MW (kDa) 30 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 200mM NaCl, 50% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAAGFKT VEPLEYRRF LKENCPRDGR ELGEFRITTV NIGSISTADG SALVKLGNTT VICGVKAFA APSTDAPDKG YVVPNVLDLP LCSSRFRSGP PGEEAQVASQ FIADVIENSQ IIQKEDLCIS PGKLWVLYC DLICLDYDGN ILDACTFALL AALKNVQLPE VTINEETALA EVNLKKKSYL NIRTHPVATS FAVFDDTLI VDPTGEEHL ATGTLTIVMD EEGKLCLHK PGGSGLTGAK LQDCMSRAVT RHKEVKKLMD EVIKSMKPK

Expression System E. coli

Purity > 90% by SDS-PAGE.


Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68251-pro Image

3μg Human EXOSC8 protein (GTX68251-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).