

Mouse CRMP1 protein, His tag


Cat. No. GTX68256-pro

Species	Mouse	Package 10 µg
Applications		
Observed MW (kDa)	62 kDa.	
Properties		
Form	Liquid	
Buffer	20mM Tris-HCl, 200mM NaCl, 40% Glycerol, 0.1mM PMSF	
Preservative	No preservatives	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Concentration	0.25 mg/ml (Please refer to the vial label for the specific concentration.)	
Region/Sequence	Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSEFMSHQG KKSIPHITSD RLLIRGGRII NDDQSFYADV YLEDGLIKQI GENLIVPGGV KTIEANGRMV IPGGIDVNTY LQKPSQGMTS ADDFFQGTKA ALAGGTTMII DHVVEPGSS LLTSFEKWHE AADTKSCCDY SLHVDITSWY DGVREELEV L VQDKGVNSFQ VYMAYKDLYQ MSDSQLYEAF TFLKGLGAVI LVHAENGDLI AQEQKRILEM GITGPEGHAL SRPEELEAEA VFRAIAIAGR INCPVYITKV MSKSAADIIA LARKKGPLVF GEPIAASLGT DGTHYWSKNW AKAAAFVTSP PLSPDPTTPD YLTSLACGD LQVTGSGHCP YSTAQKAVGK DNFTLIPEGV NGIEERMTVV WDKAVATGKM DENQFVAVTS TNAAKIFNLY PRKGRIAVGS DADVVIWDPD KMKTITAKSH KSTVEYNIFE GMECHGSPLV VISQGKIVFE DGNISVSKGM GRFIPRKFPF EHLYQVRIR SKVFGHLSVS RGMVDGPVYE VPATPKHAAP APSAKSSPSK HQPPPIRN LH QSNFSLGAQ IDDN NPRRTG HRIVAPPGGR SNITSLG	
Expression System	E. coli	
Purity	> 90% by SDS-PAGE.	
Conjugation	Unconjugated	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68256-pro Image

3μg Mouse CRMP1 protein (GTX68256-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).