

Mouse Annexin A1 protein, His tag

Cat. No. GTX68273-pro

Species	Mouse
---------	-------

Package	100 µg
---------	--------

Applications

Observed MW (kDa)	39 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	PBS, 20% Glycerol, 1mM DTT
--------	----------------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAMVSEF LKQARFLENQ EQEYVQAVKS YKGGPGSAVS PYPFNVSSD VAALHKAIMV KGVDEATIID ILTKRTNAQR QQIKAAYLQE NGKPLDEVLR KALTGHLEEV VLAMLKTPAQ FDADELRGAM KGLGTDEDTL IEILTTRSNE QIREINRVYR EELKRDIAKD ITSDTSGDFR KALLALAKGD RCQDLSVNQD LADTDARALY EAGERRKGTG VNVFTTILTS RSFPHLRVVF QNYGKYSQHD MNKALDLELK GDIEKCLTI VKCATSTPAF FAEKLYEAMK GAGTRHKALI RIMVSRSEID MNEIKVFYQK KYGISLCQAI LDETGGDYEK ILVALCGGN

Expression System	E. coli
-------------------	---------

Purity	> 95% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68273-pro Image

3μg Mouse Annexin A1 protein (GTX68273-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).