

Mouse Tyrosine Hydroxylase protein, His tag

Cat. No. GTX68296-pro

Species	Mouse
---------	-------

Package

50 µg

Applications

Observed MW (kDa)	65 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	PBS, 10% Glycerol
--------	-------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	---

Region/Sequence

Full length protein, C-terminal His-Tag; ADPMPTPSAS SPQPKGFRRA VSEQDTKQAE AVTSPRFIGR RQSLIEDARK
EREAAAAAAAA AAVASAEPPN PLEAVVFEER DGNVNLNLLF SLRGTKPSSL SRALKVFETF EAKIHHLLETR PAQRPLAGSP
HLEYFVRFEV PSGDLAALLS SVRRVSDDVR SAREDKVPWF PRKVSSELDKC HHLVTKFDPD LDLDHPGFSD QAYRQRRKL
AEIAFYKQG EPIPHVEYTK EEIATWKEVY ATLKGLYATH ACREHLEAFQ LLERYCGYRE DSIPQLEDVS HFLKERTGFQ LRPVAGLLSA
RDFLASLAFR VFQCTQYIRH ASSPMHSPEP DCHELLGHV PMLADRTFAQ FSQDIGLASL GASDEEIEKL STVWFTVEF
GLCKQNGELK AYGAGLLSSY GELLHSLSEE PEVRAFDPD AAVQPYQDQT YQPVYFVSES FSDAKDKLRN YASRIQRPFS
VKFDPYTLAI DVLDSPTIR RSLEGVQDEL HTLTQALSAI SHHHHHH

Expression System	Baculovirus
-------------------	-------------

Purity	> 90% by SDS-PAGE.
--------	--------------------

Endotoxin	< 1.0 EU/µg (determined by LAL method)
-----------	--

Conjugation	Unconjugated
-------------	--------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 04 Page 1 of 2

DATA IMAGES



GTX68296-pro Image

3μg Mouse Tyrosine Hydroxylase protein (GTX68296-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).