

# Human FUS2 protein, His tag

**Cat. No. GTX68351-pro**
**Species** Human

**Package**

10 µg

**Applications**
**Observed MW (kDa)** 31 kDa.

**Properties**
**Form** Liquid

**Buffer** 20mM MES buffer, 150mM NaCl, 30% Glycerol

**Preservative** No preservatives

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MQELTSPGP AKLTPTLDPT HRMELILSTS  
PAELTDPAC QPKLPDSTC QPEMTFNPGP TELTDPEHQ PEETPAPSLA ELTLEPVHRR PELLDACADL INDQWPRSRT  
SRLHSLGQSS DAFPLCLMLL SPHTLEAAP VVGHARLSR VLNQPSLLV ETVVVARALR GRGFGRRLME GLEVFARARG  
FRKLHLTHD QVHFYTHLGY QLGEVQGLV FTSRRLPATL LNAFPTAPSP RPPRKAPNLT AQAAPRGPKG PPLPPPPPLP ECLTISPPVP  
SGPPSKSLE TQYQNVGRGP IFWMEKDI

**Expression System** E. coli

**Purity** > 95% by SDS-PAGE.

**Conjugation** Unconjugated

**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX68351-pro Image

3μg Human FUS2 protein (GTX68351-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).