

Human BRMS1 protein, T7 tag

Cat. No. GTX68360-pro

Species Human

Package
20 µg

APPLICATION

Observed MW 28 kDa.

PROPERTIES

Form Liquid

Buffer 20mM Tris-HCl, 0.4M Urea, 10% Glycerol

Preservative No preservative

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal T7-Tag; MASMTGGQQM GRGSHMPVQP PSKDTEEMEA EGDSAAEMNG EEEEEEEERS GSQTESEES SEMDDDYER RRSECVSEML DLEKQFSELK EKLFRERLSQ LRLRLEEVEGA ERAPEYTEPL GGLQRLKIR IQVAGIYKGF CLDVIRNKYE CELQGAQHLL ESEKLLLYDT LQGELQERIQ RLEEDRQSLD LSSEWWDDKL HARGSSRSWD SLPPSKRKKA PLVSGPYIVY MLQEIDILED WTAIKKARAA VSPQKRKSDG P

Expression System E. coli

Purity > 80% by SDS-PAGE.

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68360-pro Image

3μg Human BRMS1 protein (GTX68360-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).