

Human GRASP55 protein, His tag

Cat. No. GTX68369-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 47 kDa.

Properties

Form Liquid**Buffer** PBS**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMGSSQSV EIPGGGTGEY HVLRVQENSP
GHRAGLEPFF DFIVSINGSR LNKDNDTLKD LLKANVEKPV KMLIYSSKTL ELRETSVTPS NLWGGQGGLG VSIRFCSFDG
ANENVWHVLE VESNSPAALA GLRPHSDYII GADTMVNESE DLFSLIETHE AKPLKLYVYN TDTDNCREVI ITPNSAWGGE
GSLGCGIGYG YLHRIPTRPF EEGKKISLPG QMAGTPITPL KDGFEVQLS SVNPPSLSP GTTGIEQSLT GLSISSTPPA VSSVLSTGVP
TVPLLPPQVN QSLTSVPPMN PATTLPLMP LPAGLPNLPN LNLNLPAPHI MPGVGLPELV NPGLPPLPSM PPRNLPGIAP
LPLPSEFLPS FPLVPESSA ASSGELLSSL PPTSNA PSDP ATTTAKADAA SSLTVDVTPP TAKAPTTVED RVGDSTPVSE KPVSAAVDAN
ASESP

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68369-pro Image

3μg Human GRASP55 protein (GTX68369-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).