

Human TES protein, His tag

Cat. No. GTX68373-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 48 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 100mM NaCl, 10% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MDLENKVKKM GLGHEQGFGA PCLKCKEKCE
GFELHFWRKI CRNCKCGQEE HDVLLSNEED RKGKLFEDT KYTTLIAKLK SDGIPMYKRN VMILTNPVAA KKNVSINTVT
YEWAPPVQNN ALARQYMQML PKEKQPVAGS EGAQYRKKQL AKQLPAHDQD PSKCHELSPR EVKEMEQLVK KYKSEALGVG
DVKLPCEMDA QGPKQMNIPG GDRSTPAAVG AMEDKSAEHK RTQYSCYCCK LSMKEGDP AI YAERAGYDKL WHPACFVCST
CHELLVDMYI FWKNEKLYCG RHYCDSEKPR CAGCDELIFS NEYTAENQN WHLKHFCFCD CDSILAGEIY VMVNDKPVCK
PCYVKNHAVV CQGCHNAIDP EVQRTYNNF SWHASTECFL CSCCSKCLIG QKFMPVEGMV FCSVECKKRM S

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68373-pro Image

3μg Human TES protein (GTX68373-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).