

Human SNX5 protein, His tag

Cat. No. GTX68395-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 47 kDa.

Properties

Form Liquid**Buffer** PBS, 20% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAAVPEL LQQQEEDRSK LRSVSVDLNV
DPSLQIDIPD ALSERDKVKF TVHTKTLPT FQSPEFSVTR QHEDFWLHD TLITTDYAG LIIPPAPTKP DFDGPREKMQ KLGEGECSMT
KEEFAKMKQE LEAEYLAVFK KTVSSHEVFL QRLSSHPVLS KDRNFHVLE YDQDLSVRRK NTKEMFGGFF KSVVKSADDEV
LFTGVKEVDD FFEQEKNFLI NYNRIKDSC VKADKMTRSH KNVADDYIHT AACLHSLALE EPTVIKKYLL KVAELFEKLR KVEGRVSSDE
DLKLTLLRY YMLNIEAAKD LLYRRTKALI DYENSNKALD KARLKSQDVK LAEAHQEQCC QKFEQLSESA KEELINFKRK RVAAFRKNLI
EMSELEIKHA RNNVSLQSC IDLFKNN

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx68395-pro Image

3µg Human SNX5 protein (GTx68395-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).