

Human ARFIP1 protein, His tag

Cat. No. GTX68400-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 42 kDa.

Properties

Form Liquid**Buffer** PBS, 30% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMASQSPK NSAAEIPVTS NGEVDDSDREH
SFNRDLKHSL PSGGLSETQ ITSHGFDNTK EGVIEAGAFQ GGQRTQTKSG PVILADEIKN PAMEKLELVR KWSLNTYKCT
RQIIEKLGR GSRTVDLELE AQIDILRDNK KKYENILKLA QTLSTQLFQM VHTQRQLGDA FADLSLSLE LHEEFGYNAD
TQKLLAKNGE TLLGAINFFI ASVNTLVNKT IEDTLMVKQ YESARIEYDA YRTDLEELNL GPRDANTLPK IEQSQHLEFQA
HKEKYDKMRN DVSVKLFLE ENKVVLHNQ LVLFHNAIAA YFAGNQKQLE QTLKQFHIKL KTPGVDAPSW LEEQ

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68400-pro Image

3μg Human ARFIP1 protein (GTX68400-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).