

# Human MAT2B protein

**Cat. No. GTX68416-pro**

**Species** Human

**Package**  
100 µg

## Applications

**Observed MW (kDa)** 38 kDa.

## Properties

**Form** Liquid

**Buffer** 20mM Tris-HCl, 10% Glycerol, 1mM DTT

**Preservative** No preservatives

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)

**Region/Sequence** Full length protein, MPMPEDMEQ EEVNIPNRRV LVTGATGLLG RAVHKEFQQN NWHAVGCGFR RARPKFEQVN LLDNSNAVHHI IHDFQPHVIV HCAAERRPDV VENQPDAAASQ LNVDAAGNLA KEAAAVGAFL IYSSDYVFD GTNPPYREED IPAPLNLYGK TKLDGEKAVL ENNLGA AVL IPILYGEVEK LEESAVTVMF DKVQFSNKSA NMDHWQQRFP THVKDVATVC RQLAEKRMLD PSIKGTFHWS GNEQMTKYEM ACAIADAFNL PSSHLRPITD SPVLGAQRPR NTQLDCSKLE TLGIGQRTPF RIGIKESLWP FLIDKRWRQT VFH

**Expression System** E. coli

**Purity** > 90% by SDS-PAGE.

**Conjugation** Unconjugated


**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX68416-pro Image

3μg Human MAT2B protein (GTX68416-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).