

Human CRNN protein, His tag

Cat. No. GTX68456-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 54 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSSLVPRGSH MPQLQNING IIEAFRRYAR TEGNCTALTR
GELKRLLQE FADVIVKPHD PATVDEVRLR LDEDHTGTVE FKEFLVLVFK VAQACFKTLS ESAEGACGSQ ESGSLHSGAS
QELGEGQSRG TEVGRAGKGQ HYGSSSHRQS QQGSRGQNRG GVQTQGGATG SAWVSSYDRQ AESQSQERIS PQQLSGQTE
QTQKAGEGKR NQTTEMRPER QPQTRQDRA HQTGETVTGS GTQTQAGATQ TVEQDSSHQT GRTSKQTQEA TNDQNRGTET
HGQGRSQTSG AVTGGHAQIQ AGHTQTPTQ TVEQDSSHQT GSTSTQTQES TNGQNRGTEI HGQGRSQTSG AVTGGHTQIQ
AGSHTETVEQ DRSQTVSHGG AREQGQTQTQ PGSGQRWMQV SNPEAGETVP GGQAQTGAST EPGRQEWSSST HPRRCVTEGQ
GDRQPTVVGE EWVDDHSRET VILRLDQGNL HTSVSSAQGQ DAAQSEEKRG ITARELYSYL RSTKP

Expression System E. coli

Purity > 85% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES



GTX68456-pro Image

3μg Human CRNN protein (GTX68456-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).