

Human EXOSC3 protein, His tag

Cat. No. GTX68463-pro

Species Human

Package
10 µg

Applications

Observed MW (kDa) 30 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 150mM NaCl, 30% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAEPASV AAESLAGSRA RAARTVLGQV VLPGEELLPL EQEDAEGPGG AVERPLSLNA RACSRVRVC GPGLRRCGDR LLVTKCGRRLR HKEPGSGSGG GVVWVDSQQK RYVPVKGDHV IGIVTAKSGD IFKVDVGGSE PASLSYLSFE GATKRNRPNV QVGDLIYGQF VWANKDMEPE MVCIDSCGRA NGMGVIGQDG LLFKVTLGLI RLLAPDCEI IQEVGKLYPL EIVFGMNGRI WVKAKTIQQT LILANILEAC EHMTSDQRKQ IFSRLAES

Expression System E. coli

Purity > 85% by SDS-PAGE.

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68463-pro Image

3μg Human EXOSC3 protein (GTX68463-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).