

Human SLA/LP protein, His tag

Cat. No. GTX68473-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 56 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 150mM NaCl, 40% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MNRESFAAGE RLVSPAYVRQ GCEARRSHEH
LIRLLLEKGGK CPENGWDEST LELFLHELAI MDSNNFLGNC GVGEREGRVA SALVARRHYR FIHGIGRSGD ISAVQPKAAG
SSLLNKITNS LVLDIKLAG VHTVANCFFV PMATGMSLT L CFTLRHKRP KAKYIWPRI DQKSCFKSMI TAGFEPWIE NVLEGDELRT
DLKAVEAKVQ ELGPDCLICI HSTTSCFAPR VPDRLLELAV ICANYDIPHI VNNAYGVQSS KCMHLIQQGA RVGRIDAFVQ
SLDKNFMVPV GGAIAGFND SFIQEISKMY PGRASASPSL DVLITLLSLG SNGYKLLKE RKEMFSYLSN QIKKLSEAYN ERLHPTHNP
ISLAMTLKTL DEHRDKAVTQ LGSMLFTRQV SGARWPLGS MQTVSGYTFR GFMSHTNNYP CAYLNAASAI GMKMQDQDLF
IKRLDRCLKA VRKERSKESD DNYDKTEDVD IEEMALKLDN VLLDTYQDAS S

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68473-pro Image

3µg Human SLA/LP protein (GTX68473-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).