

Human VRK3 protein, His tag

Cat. No. GTX68490-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 53 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 150mM NaCl, 40% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMISFCPD CGKSIQAAFK FCPYCGNSLP VEEHVGSTQF VNPVVSFQG SKRGLNSSFE TSPKKVKWSS TVTSPRLSLF SDGDSESED TLSSERSKG SGRPPTPKS SPQKTRKSPQ VTRGSPQKTS CSPQKTRQSP QTLKRSRVTT SLEALPTGTV LTDKSGRQWK LKSFQTRDNQ GILYEAAPTS TLTCDSGPQK QKFSCLKDAK DGRLFNEQNF FQRAAKPLQV NKWKLYSTP LLAIPTCMGF GVHQDKYRFL VLPSLGRSLQ SALDVSPKHV LSERSVLQVA CRLDALEFL HENEYVHGNV TAENIFVDPE DQSQVTLAGY GFAFRYCPSG KHVAYVEGSR SPHEGDLEFI SMDLHKGCGP SRRSDLQSLG YCMLKWLYGF LPWTNCLPNT EDIMKQKQKF VDKPGPFVGP CGHWIRPSET LQKYLKVVMA LTYEEKPPYA MLRNNLEALL QDLRVSPYDP IGLPMVP

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68490-pro Image

3μg Human VRK3 protein (GTX68490-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).