

Human NANS protein, His tag

Cat. No. GTX68545-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 40 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 10% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MPLELELCPG RWVGGQHPCF IIAEIGQNHQ
GDLDAKRMIRMAKECGADC AKFQKSELEF KFNKALERP YTSKHSWGKT YGEHKKHLEF SHDQYRELQR YAEVGIFFT
ASGMDEMAVE FLHNLNVPFF KVGSGDTNNF PYLEKTAKKG RPMVISSGMQ SMDTMKQVYQ IVKPLNPNFC FLQCTSAAYPL
QPEDVNLRLVI SEYQKLFPDI PIGYSGHETG IAISVAAVAL GAKVLERHIT LDKTWKGS DH SASLEPGELA ELVRSVRLVE RALGSPTKQL
LPCMACNEK LGKSVVAKVK IPEGTILTMD MLTVKVGEPK GYPEDIFNL VGKKVLVTVE EDDTIMEELV DNHGKKIKS

Expression System E. coli

Purity > 95% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES



GTX68545-pro Image

3μg Human NANS protein (GTX68545-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).