

Human DDX56 protein, His tag

Cat. No. GTX68554-pro**Species** Human**References** (1)**Package**

100 µg

Applications**Observed MW (kDa)** 62 kDa.**Properties****Form** Liquid**Buffer** 20mM Tris-HCl, 0.4M Urea, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMEDSEAL GFEHMGLDPR LLQAVTDLGW
 SRPTLIQEKA IPLALEGKDL LARARTGSGK TAAYAIPMLQ LLLHRKATGP VVEQAVRGLV LVPTKELARQ AQSMIQQLAT
 YCARDVRVAN VSAAEDSVSQ RAVLMEKPDV VVGTPSRILS HLQQDSLKLR DSLELLWDE ADLLFSFGFE EELKSLLCHL
 PRIYQAFILMS ATFNEDVQAL KELILHNPVT LKLQESQLPG PDQLQQFQV CETEEKFLL LYALLKLSLI RGKSLLFVNT LERSYRLRLF
 LEQFSIPTCV LNGELPLRSR CHIISQFNQG FYDCVIATDA EVLGAPVKGK RRGRGPKGDK ASDPEAGVAR GIDFHHVSAV
 LNFDLPPTPE AYIHRAGRTA RANNPVILT FVLPTQFHL GKIEELLSGE NRGPIILPYQ FRMEEIEGFR YRCRDAMRSV TKQAIREARL
 KEIKEELLHS EKLKTYFEDN PRDLQLLRHD LPLHPAVWKP HLGHVDPDYL PPALRGLVRP HKKRKKLSSS CRKAKRAKSQ
 NPLRSFKHKG KKFRPTAKPS

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and
 publications, please visit our [website](#).

DATA IMAGES



GTx68554-pro Image

3μg Human DDX56 protein (GTx68554-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).