

Human FBLIM1 protein, His tag

Cat. No. GTX68558-pro

Species Human

Package

20 µg

Applications

Observed MW (kDa) 41 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMASKPEK RVASSVFITL APPRRDVAVA
EEVRQAVCEA RRGPRWEAPA PMKTPEAGLA GRPSPWTTTPG RAAATVPAAP MQLFNGGCPP PPPVLDGEDV LPDLDLLPPP
PPPPVLLPS EEEAPAPMGA SLIADLEQLH LSPPPPPQA PAEGPSVQPG PLRPMEEELP PPPAEPVEKG ASTDICAFCF KTVSPRELAV
EAMKRQYHAQ CFTCRTRRRQ LAGQSFYQKD GRPLCEPCYQ DTLERCCKGCG EVVRDHIIRA LGQAFHPSCF TCVTCARCIG
DESFALGSQN EVYCLDDFYR KFAPVCSICE NPIIPRDGKD AFKIECMGRN FHENCYRCED CRILLSVEPT DQGCYPLNNH
LFCKPCHVKR SAAGCC

Expression System E. coli

Purity > 85% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68558-pro Image

3μg Human FBLIM1 protein (GTX68558-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).