

Human PGM2 protein, His tag

Cat. No. GTX68588-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 68 kDa.

Properties

Form Liquid**Buffer** PBS, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAAPEGS GLGEDARLDQ ETAQWLRWDK
NSLTLEAVKR LIAEGNKEEL RLCFGARMEF GTAGLRAAMG PGISRMNDLT IIQTTQGFCR YLEKQFSDLK QKGIVISFDA
RAHPSSGGSS RRFARLAATT FISQGIPVYL FSDITPTPFV PFTVSHLKC AGIMITASHN PKQDNGYKVY WDNGAQIISP HDKGISQAIE
ENLEPWPQAW DDSLIDSSPL LHNPSASINN DYFEDLKCYC FHRSVNRETK VKFVHTSVHG VGHSEVQSAF KAFDLVPPEA
VPEQKDPDPE FPTVKYPNPE EGKGVLTLSF ALADTKARI VLANDPDADR LAVAQKQDSG EWRVFSGNEL GALLGWWLFT
SWKEKNQDRS ALKDTYMLSS TVSSKILRAI ALKEGFHFEE TLGFKWMGN RAKQLIDQGK TVLFAFEEAI GYMCCPFVLD
KDGVSAAVIS AELASFLATK NLSLSQQLKA IYVEYGYHIT KASYFICHQD ETIKKLFENL RNYDGKNNYP KACGKFEISA IRDLTTGYDD
SQPDKKAVLP TSKSSQMITF TFANGGVATM RTSGTEPKIK YYAELCAPPG NSDPEQLKKE LNELVSAIEE HFFQPQKYNL QPKAD

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68588-pro Image

3μg Human PGM2 protein (GTX68588-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).