

Human TBC1D22B protein, His tag

Cat. No. GTX68601-pro

Species Human

Package

10 µg

Applications

Observed MW (kDa) 59 kDa.

Properties

Form Liquid

Buffer PBS, 30% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAAENSK QFWKRSALP GSIQPVYGAQ HPPLDPRLTK NFIKERSKVN TVPLKNKKAS SFHEFARNTS DAWDIGDDEE EDFSSPSFQT LNSKVALATA AQVLENHSL RVKERSQST TSDVPANYKV IKSSDAQLS RNSSDTCLRN PLHKQQLPL RPIIPLVARI SDQNASGAPP MTVREKTRLE KFRQLSSQN TDLDELKCS WPGVPREVRP ITWRLLSGYL PANTERRKLT LQRKREEYFG FIEQYYDSRN EEHHQDITYRQ IHIDIPRTNP LIPLFQQPLV QEIFERILFI WAIRHPASGY VQGINDLVTP FFVFLSEYV EEDVENFDVT NLSQDMLRSI EADSFWCMSK LLDGIQDNYT FAQPGIQKKV KALEELVSRI DEQVHNHFR YEVEYLQFAF RWMNNLLMRE LPLRCTIRLW DTYQSEPEGF SHFHLVCAAA FLIKWRKEIL DEEDFQGLLM LLQNLPTIHW GNEEIGLLA EAYRLKYMFA DAPNHYRR

Expression System E. coli

Purity > 80% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES



GTX68601-pro Image

3μg Human TBC1D22B protein (GTX68601-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).