

Human PBK protein, His tag

Cat. No. GTX68621-pro
Species Human

Package
100 µg

Applications
Observed MW (kDa) 36 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 10% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSHMEGISH FKTPSKLSEK KKSVCSTPT INIPASPFMQ KLGFGTGVNV YLMKRSPRGL SHSPWAVKKI NPICNDHYRS VYQKRLMDEA KILKSLHHPN IVGYRAFTEA NDGSLCLAME YGGEKSLNDL IEERYKASQD PFPAAILKV ALNMARGLKY LHQEKLLHG DIKSSNVVIK GDFETIKICD VGVSLPLDEN MTVTDPEACY IGTEPWKPKE AVEENGVID KADIFAFGLT LWEMMTLSIP HINLSNDDDD EDKTFDESDF DDEAYYAALG TRPPINMEEL DESYQKVIEL FSVCTNEDPK DRPSAAHIVE ALETDV

Expression System E. coli

Purity > 95% by SDS-PAGE.


Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68621-pro Image

3μg Human PBK protein (GTX68621-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).