

Human THAP11 protein, His tag

Cat. No. GTX68653-pro

Species	Human
---------	-------

Package

50 µg

Applications

Observed MW (kDa)	34 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 0.4M Urea, 10% Glycerol
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMPGFTCC VPGCYNNNSHR DKALHFYTFP
KDAELRRLWL KNVSRAGVSG CFSTFQPTTG HRLCSVHFQG GRKTYTVRVP TIFPLRGVNE RKVARRPAGA AAARRRQQQQ
QQQQQQQQQQ QQQQQQQQQQ QQQQQSSPSA STAQTAQLQP NLVSASAAVL LTLQATVDSS QAPGSVQPAP ITPTGEDVKP
IDLTVQVEFA AAEGAAAAAA ASELQAATAG LEAAECPMGP QLVVVGEEGF PDTGSDHSYS LSSGTTEEEL LRKLNEQRDI
LALMEVKMKE MKGSIRHLRL TEAKLREELR EKDRLLAMAV IRKKHGM

Expression System	E. coli
-------------------	---------

Purity	> 80% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES



GTX68653-pro Image

3μg Human THAP11 protein (GTX68653-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).