

Human NDRG2 protein, His tag

Cat. No. GTX68661-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 41 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 100mM NaCl, 400mM Urea, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSHMAELQE VQITEEKPLL PGQTPEAAKT
HSVETPYGSV TFTVYGTPKP KRPAILTYHD VGLNYKSCFQ PLFQFEDMQE IIQNFVRVHV DAPGMEEGAP VFPLGYQYPS
LDQLADMIPC VLQYLNFSIT IGVGVGAGAY ILARYALNHP DTVEGLVLIN IDPNAKGWMD WAAHKLTGLT SSIPEMILGH
LFSQEELSGN SELIQYRNII ITHAPNLDNI ELYWNSYNNR RDLNFERGGD ITRCPVMLV VGDQAPHEDA VVECNSKLDP
TQTSFLKMAD SGGQPQLTQP GKLTEAFKYF LQGMGYMASS CMTRLRSRT ASLTAASVD GNRSRRTLS QSSESGTLSS
GPPGHTMEVS C

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68661-pro Image

3μg Human NDRG2 protein (GTX68661-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).