

Human GBA3 protein, His tag

Cat. No. GTX68664-pro

Species	Human
---------	-------

Package	100 µg
---------	--------

Applications

Observed MW (kDa)	54 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	PBS, 20% Glycerol, 1mM DTT
--------	----------------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	---

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAPAGF GWAAATAAYQ VEGGWDADGK GPCVWDTFTH QGGERVFNQ TGDVACGSYT LWEEDLKCIK QLGLTHYRFS LWSRLLPDG TTGFNQKI DYYNKIIDDL LKNGVTPVT LYHFDLPQTL EDQGGWLSEA IIESFDKYAQ FCFSTFGDRV KQWITINEAN VLSVMSYDLG MFPPGIPHFG TGGYQAAHNL IKAHARSWHS YDSLFRKKQK GMVSLSLFAV WLEPADPNSV SDQEAARKAI TFHLDLFAKP IFIDGDYPEV VKSQIASMSQ KQGYPSSRLP EFTEEEKMI KGTADFFAVQ YYTTRLIKYQ ENKKGELGIL QDAEIEFFPD PSWKNVDWIY VVPWGVCKLL KYIKDTYNNP VIYITENGFP QSDPAPLDDT QRWEYFRQTF QELFKAIQLD KVNQVYCAW SLLDNFEWNQ GYSSRFLFH VDFEDPARPR VPYSAKEYA KIIRNNGLEA HL

Expression System	E. coli
-------------------	---------

Purity	> 85% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX68664-pro Image

3μg Human GBA3 protein (GTX68664-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).