

Mouse Visfatin protein, His tag


Cat. No. GTX68676-pro

Species	Mouse	Package 100 µg
Applications		
Observed MW (kDa)	55 kDa.	
Properties		
Form	Liquid	
Buffer	PBS	
Preservative	No preservatives	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)	
Region/Sequence	Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MNAAAEAEFN ILLATDSYKV THYKQYPPNT SKVSYFECEK EKTENSKVR KVKYEETVFY GLQYILNKYL KGKVVTKEDI QEAKEVYREH FQDDVFNERG WNYILEKYDG HLPYEVKAVP EGSVIPRGNV LFTVENTDPE CYWLTNWIET ILVQSWYPIT VATNSREQKK ILAKYLLETS GNLDGLEYKL HDFGYRGVSS QETAGIGASA HLVNFKGTD T VAGIALIKKY YGKDPVPGY SVPAAEHSTI TAWGKDHEKD AFEHIVTQFS SVPVSVWSDS YDIYNACEKI WGEDLRHLIV SRSTEAPLII RPDSGNPLDT VLKVL DILGK KFPVTENSKG YKLLPPYLRV IQGDGVDINT LQEIVEGMKQ KKWSIENVSF GSGGALLQKL TRDLLNCSFK CSYVVTNGLG VNVFKDPVAD PNKRSKKGRL SLHRTAGNF VTLEEGKGD L EEYGHDL LHT VFKNGKVTKS YSFDEVKNA QLNIEQDVAP H	
Expression System	E. coli	
Purity	> 90% by SDS-PAGE.	
Endotoxin	< 1.0 EU/µg (determined by LAL method)	
Conjugation	Unconjugated	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68676-pro Image

3μg Mouse Visfatin protein (GTX68676-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).