

Human CPEB1 protein, His tag

Cat. No. GTX68699-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 63 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 0.4M Urea, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAPLEE EAGRIKDCWD NQEAPALSTC
SNANIFRRIN AILDNSLDFS RVCTTPINRG IHDHLPDFQD SEETVTSRML FPTS AQESSR GLPDANDLCL GLQSLSLTGW
DRPWSTQDSD SSAQSSTHSV LSMLHNPLGN VLGKPLSFL PLDPLGSDLV DKFPAPSVRG SRLDTRPILD SRSSPSDSD
TSGFSSGSDH LSDLISSLRI SPPLPFLSLS GGGPRDPLKM GVGS RMDQEQ AALAAVTPSP TSASKRWPGA SVWPSWDLLE
APKDPFSIER EARLHRQAAA VNEATCTWSG QLPPRNYKNP IYSCKVFLGG VPWDITEAGL VNTFRVFGSL SVEWPGKD GK
HPRCPPKGYV YLVFELEKSV RSLAQACSHD PLSPDGLSEY YFKMSSRRMR CKEVQVIPWV LADSNFVRSP SQRLDPSRTV
FVGALHGMLN AEALAILND LFGGVYAGI DTDKHKYPIG SGRVTFNNQR SYLKAVSAAF VEIKTKFTK KVQIDPYLED
SLCHICSSQP GPFFCRDQVC FKYFCRSCWH WRHSM EGLRH HSPLMRNQKN RDSS

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68699-pro Image

3μg Human CPEB1 protein (GTX68699-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).