

Mouse Elongin B protein, His tag

Cat. No. GTX68712-pro
Species Mouse

Package

50 µg

Applications
Observed MW (kDa) 13 kDa.

Properties
Form Liquid

Buffer PBS, 1mM DTT, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMDVFLMI RRHKTTIFTD AKESSTVFEL KRIVEGILKR PPEEQRLYKD DQLDDGKTL GECGFTSQTA RPQAPATVGL AFRADDTFEA LRIEPFSSPP ELPDVMKPQD SGGANEQAV Q

Expression System E. coli


Purity > 95% by SDS-PAGE.

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX68712-pro Image

3µg Mouse Elongin B protein (GTX68712-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.


For full product information, images and publications, please visit our [website](#).