

Human TUG protein, His tag

Cat. No. GTX68725-pro
Species Human

Package

20 µg

Applications
Observed MW (kDa) 60 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAAPAGG GGSVSVLAP NGRRHVTKVT
PSTVLLQVLE DTCRRQDFNP CEYDLKFQRS VLDLSLQWRF ANLPNNAKLE MVPASRSREG PENMVRIALQ LDDGSRLQDS
FCSGQTLWEL LSHFPQIREC LQHPGGATPV CVYTRDEV TG EAALRGTTLQ SLGLTGGSAT IRFVMKCYDP VGKTPGSLGS
SASAGQAAAS APLPLESGEL SRGDLSPED ADTSGPCCEH TQEKQSTRAP AAAPFVPFSG GGQRLGGPPG PTRPLTSSSA
KLPKSLSPG GPSKPKSKS GQDPQQEQEQ ERERDPQQEQ ERERPVDREP VDREPVVCHP DLEERLQAWP AELPDEFFEL
TVDDVRRRLA QLKSERKRL EAPLVTKAFR EAQIKEKLER YPKVALRVLF PDRYVLQGFF RPSETVGLR DFVRSHLGNP ELSFYLFITP
PKTVLDDHTQ TLFQANLFPA ALVHLGAEAP AGVYLEPGLL EHAISPSAAD VLARYMSRA AGSPSPLPAP DPAPKSEPAA
EEGALVPPPEP IPGTAQPVKR SLGKVPKWLK LPASKR

Expression System E. coli

Purity > 85% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68725-pro Image

3μg Human TUG protein (GTX68725-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).