

Human PA1 protein, His tag

Cat. No. GTX68732-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 28 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 10% Glycerol, 2mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MSLARGHGD T AASTAAPLSE EGEVTSG LQA LAVEDTGGPS ASAGKAED EG GREETER E GSGGEEAQGE VPSAGGEEPA EEDSEDWCVP CSDEEVELPA DGQPWMPPPS EIQRLYELLA AHGTLELQAE ILPRRPPTPE AQSEEEERSDE EPEAKEEEEE KPHMPTEFDF DDEPVTPKDS LIDRRRTPGS SARSQKREAR LDKVLSDMKR HKKLEEQILR TGRDLFSLDS EDPSPASPPL RSSGSSLFPR QRKY

Expression System E. coli

Purity > 85% by SDS-PAGE.

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68732-pro Image

3μg Human PA1 protein (GTX68732-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).