

Human PPCS protein, His tag

Cat. No. GTX68742-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 34 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MAEMDPVAEF QPPGAARWA EVMARFAARL GAQGRRVVLV TSGGTKVPLE ARPVRFLDNF SSGRRGATSA EAFLAAGYGV LFLYRARSF PYAHRFPPT WLSALRPSGP ALSGLLSLEA EENALPGFAE ALRSYQAAAA AGTFLAVEFT TLADYLHLLQ AAAQALNPLG PSAMFYLAAS VSDFYVPVSE MPEHKIQSSG GPLQITMKMV PKLLSPLVKD WAPKAFIISF KLETDPALVI NRARKALEIY QHQVVANIL ESRQSFVIV TKDSETKLL SEEEIEKGVE IEEKIVDNLQ SRHTAFIGDR N

Expression System E. coli

Purity > 95% by SDS-PAGE.

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68742-pro Image

3μg Human PPCS protein (GTX68742-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).