

Hec1 (non-phospho Ser 165) blocking peptide

Cat. No. GTX70012-PEP

Applications	Neutralizing /Inhibition
Species	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Neutralizing /Inhibition	Assay dependent
--------------------------	-----------------

Not tested in other applications.

Properties

Form	Liquid
-------------	--------

Buffer	Water
---------------	-------

Preservative	No preservatives
---------------------	------------------

Storage	Store at -20°C. Avoid repeated freeze/thaw cycles.
----------------	--

Concentration	2.5 mg/ml (Please refer to the vial label for the specific concentration.)
----------------------	--

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).