

APE1 antibody

Cat. No. GTX70129

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

For WB: Use at a dilution of 1:2,500. Predicted molecular weight: 38 kDa. Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW 36 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Full length human APE1 purified from E. coli
Purification	Unpurified From polyclonal serum
Conjugation	Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX70129 WB Image

GTX70129 at a 1:2500 dilution staining APE / REF 1 in A) HeLa whole cell lysate and B) Raji whole cell lysate in Western Blot (ECL).



For full product information, images and publications, please visit our [website](https://www.genetex.com).