

DNA polymerase delta antibody

Cat. No. GTX70151

Host	Mouse
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

References (1)

Package

50 µl

Applications

Application Note

For WB: Use at a dilution of 1:500. Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW 124 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Unpurified Serum
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant protein encoding amino acids 1-140 of the large subunit of human DNA polymerase d.
Purification	Unpurified From polyclonal serum.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX70151 WB Image

Detection of human polymerase delta protein in Jurkat whole cell lysate using GeneTex DNA polymerase delta polyclonal antiserum GTX70151.



For full product information, images and publications, please visit our [website](#).