

E2F1 antibody [18E10]

Cat. No. GTX70165

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse

References (2)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IP	Assay dependent

Not tested in other applications.

Calculated MW 47 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The region encoding the DNA binding domain (amino acids 244-264) of E2F-1 expressed in E. coli.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTx70165 IHC-P Image

E2F1 antibody [18E10] detects E2F1 protein at nucleus by immunohistochemical analysis.

Sample: Paraffin-embedded mouse ovary.

E2F1 stained by E2F1 antibody [18E10] (GTx70165) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx70165 WB Image

Sample (30 ug of whole cell lysate)

A: 293T

10% SDS PAGE

GTx70165 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).