

UBE3A antibody

Cat. No. GTX70169

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC
Reactivity	Human

Package

50 µl

Applications

Application Note

Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW 101 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Unpurified Serum
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant protein encoding amino acids residues 1-851 of human UBE3A expressed in E. coli.
Purification	Rabbit Polyclonal Antiserum Antiserum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX70169 IHC Image



For full product information, images and publications, please visit our [website](#).