

FEN1 antibody

Cat. No. GTX70184

Host	Mouse
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µl

Applications

Application Note

For WB: Use at an assay dependent dilution. Predicted molecular weight: 43 kDa. Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW 43 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant full-length human FEN-1 protein.
Purification	Unpurified Antiserum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX70184 WB Image

Detection of human FEN1 protein in HeLa whole cell lysate using GeneTex polyclonal antiserum GTX70184



For full product information, images and publications, please visit our [website](https://www.genetex.com).