

SAFB antibody [6F7]

Cat. No. GTX70197

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IP, IHC
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/ml
ICC/IF	Assay dependent
IP	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 103 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to amino acids 345-357 (PEARDSKEDGRKF) of the human SAFB protein.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTXT0197 WB Image


HET/SAF-B antibody [6F7] detects HET/SAF-B protein by western blot analysis.

- A. 30 µg 293T whole cell lysate/extract
- B. 30 µg A431 whole cell lysate/extract
- C. 30 µg HeLa whole cell lysate/extract
- D. 30 µg HepG2 whole cell lysate/extract

5% SDS-PAGE


HET/SAF-B antibody [6F7] (GTXT0197) dilution: 1:1000

The HRP-conjugated anti-mouse IgG antibody (GTXT21311-01) was used to detect the primary antibody.



GTXT0197 WB Image

Detection of human HET SAF-B protein in HeLa whole cell lysate using GeneTex monoclonal antibody GTXT0197



GTXT0197 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using HET/SAF-B [6F7](GTXT0197) antibody at 1:200 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).