

Jab1 antibody [8H8.5]

Cat. No. GTX70207

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IP
Reactivity	Human, Mouse, Drosophila

References (8)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IP	Assay dependent

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	The full-length mouse JAB-1 gene expressed in E. coli.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX70207 WB Image

Detection of human Jab1 protein in HeLa whole cell lysate using GeneTex monoclonal antibody GTX70207


GTX70207 WB Image

Part of a customer review using GTX70207


GTX70207 WB Image

Customer Review using GTX70207 - D10 T cells were lysed in IP buffer (20mM Tris-HCl, pH 7.5, 50mM NaCl, 50mM NaF, 10mM -glycerophosphate, 1% Triton X-100, 1mM PMSF).



For full product information, images and publications, please visit our [website](https://www.genetex.com).