

Nuclear Matrix Protein p84 antibody [5E10] (HRP)

Cat. No. GTX70220-01

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB
Reactivity	Human, Mouse

References (157)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000

Not tested in other applications.

Calculated MW 76 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.05M Tris, 0.15M NaCl, 1%BSA
Preservative	0.005% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.3 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Amino acids 15-374 of human p84 expressed in E. coli.
Purification	Affinity purified by Protein G.
Conjugation	Horseradish peroxidase(HRP)

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX70220-01 WB Image


p84 antibody [5E10] (HRP) detects p84 protein by western blot analysis.

A. 30 μ g HeLa whole cell lysate/extract

B. 30 μ g HeLa nuclear lysate/extract

7.5 % SDS-PAGE

p84 antibody [5E10] (HRP) (GTX670220-01) dilution: 1:1000



GTX70220-01 WB Image

p84 antibody [5E10] (HRP) detects p84 protein by western blot analysis.

A. 30 μ g Jurkat whole cell lysate/extract

B. 30 μ g Raji whole cell lysate/extract

C. 30 μ g NCI-H929 whole cell lysate/extract

7.5 % SDS-PAGE

p84 antibody [5E10] (HRP) (GTX70220-01) dilution: 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).