

TRIP230 antibody

Cat. No. GTX70252

Host	Mouse
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human

Package

50 µl

Applications

Application Note

For IP: Use at an assay dependent dilution. For WB: Use at a dilution of 1:5,000. Predicted molecular weight: 218 kDa. Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW 228 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum diluted with PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Amino acids 1099-1372 of Trip230 expressed in E. coli.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX70252 WB Image

Various whole cell extracts (30 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with TRIP230 antibody (GTX70252) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).