

FEN1 antibody

Cat. No. GTX70257

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

For WB: Use at a dilution of 1:2,500. Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW 43 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum diluted with PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant full-length human FEN-1 protein
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX70257 WB Image

Western blot analysis of FEN1 expression in Raji (A) and NIH-3T3 (B) whole cell lysates.



For full product information, images and publications, please visit our [website](#).