

PCNA antibody [339]

Cat. No. GTX70298

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:500
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 29 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant PCNA interdomain connector loop protein
Purification	Protein G purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX70298 WB Image

Sample (30 ug of whole cell lysate)


A: 293T

B: NIH-3T3


C: PC-12

12% SDS PAGE

GTX70298 diluted at 1:1000



GTX70298 WB Image



GTX70298 WB Image

Sample (whole cell lysate)

A: 293T 20ug

B: 293T 10ug

C: 293T 5ug

12% SDS PAGE

GTX70298 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).