

SFRP5 antibody

Cat. No. GTX70353

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µg

Applications

Application Note

For WB: Use at a dilution of 1:500 - 1:1,000 dilution. Optimal dilutions/concentrations should be determined by researcher.

Calculated MW	36 kDa. (Note)
Product Note	Residues 25-38 [APARCEEYDYYGWQ] of Human SFRP5

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Residues 25-38 [APARCEEYDYYGWQ] of human SFRP5
Purification	Affinity purified IgG
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx70353 WB Image

Western blot analysis of human SFRP5 in pancreas lysate.



For full product information, images and publications, please visit our [website](#).