

Acinus antibody

Cat. No. GTX70370

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µg

Applications

Application Note

For WB: Use at a dilution on 1:500-1:1,000. Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW	152 kDa. (Note)
Product Note	C terminus of Acinus.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Immunogen sequence information is commercially sensitive and details are not available.
Purification	Affinity purified IgG
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX70370 WB Image

Western blot analysis of Acinus in transfected COS-7 cell lysate.



For full product information, images and publications, please visit our [website](#).